SP FTS Ultra-Low

Wide Temperature Range Recirculating Chillers

Features & Benefits

- Compact chiller designed for continuous operation
- Temperature range of -90 °C to +150 °C (with auxillary heater)
- Fast ramping of process and fast cool down rates
- Optimal heat removal capacity with quick ramping capability to decrease test time
- Mechanical refrigeration meets requirements for a large variety of applications without need for expendable cryogens such as liquid nitrogen or dry ice
- Precise temperature control and stability
- RS232 interface allows remote control and data collection
- Portable, caster mounted
- Centrifugal or positive displacement pump options available




SPFTS

Fluid Selection Tip

For best pumping or stirring results, choose a fluid medium with a kinematic viscosity of 20 centistokes or less over the full operating temperature range. Remember, the heat exchanger in the recirculating cooler may be several degrees colder than your process.



SP FTS Ultra-Low



Additional Options & Accessories

High Pressure Pump

• Internal high pressure pump option available on all models can deliver 4 gpm up to 8 psi

External High Pressure Bypass

• The External High Pressure Bypass allows the user to limit maximum process pressure; mounted externally between the chiller inlet and outlet; includes 0-160 psi gauge (minimum setting 5 psi)

Digital Deluxe Controller

 Provides additional safeties, including Low Fluid Level, low Flow and Over-Temp with contact closure for external alarms; also includes external RTD for remote temperature control and monitoring, and auxiliary heater for rapid transition to +75 °C

Model		RC210	RC211	RC311	RC311GWT
Heat Removal (watts/BTUs)	-20°C -30°C -40°C -50°C -60°C -70°C -80°C	1350/4610 1250/4260 1150/3925 1000/3410 880/3005 550/1877	1350/4610 1250/4260 1150/3925 1000/3410 880/3005 550/1877	3400/11590 3000/10230 2700/9207 2400/8184 1800/6130 1200/4090 500/1705	3400/11590 3000/10230 2700/9207 2400/8184 1800/6130 1200/4090
Temperature Range °C		-80 to -10 °C	-80 to 75 °C	-90 to 75 °C	-75 °C to 150 °C
Reservoir	Liters	5	5	5	11.6
Compressor		2 @ 1hp	2 @ 1hp	2 @ 4hp	2 @ 4hp
Controller Interface (RS232)		Digital / RS232	Deluxe Digital / RS232	Deluxe Digital / RS232	Deluxe Digital / RS232
Control Stability ¹	±°C	0.1	0.1	0.3	0.3
Heater (Nominal 240V)		-	1500W	1500W	3500W
Pump: Standard - Centrifugal Option: Positive Displacement		2 to 3 gpm @ 22 psi 4 gpm up to 85 psi	2 to 3 gpm @ 22 psi 4 gpm up to 85 psi	2 to 3 gpm @ 22 psi 4 gpm up to 85 psi	2 to 3 gpm @ 22 psi 4 gpm up to 85 psi
Dimensions (w x d x h)	cm inches	53 x 74 x 122 21 x 29 x 48	53 x 74 x 122 21 x 29 x 48	53 x 74 x 145 21 x 29 x 57	53 x 74 x 145 21 x 29 x 57
Weight	kg / lbs	177 / 390	177 / 390	227 / 500	227 / 500
Electrical Configuration	60Hz Option 50 Hz Option ²	208 or 230 V 220 V	208 or 230 V 220 V	208 or 230 V, 3 Φ 220 V, 3 Φ	208 or 230 V, 3 Φ -
Utility Requirements	Fluid Connections Water Connections ³	1/2-in NPT, male 1/2-in NPT, female	1/2-in NPT, male 1/2-in NPT, female	1/2-in NPT, male 1/2-in NPT, female	1/2-in NPT, male 1/2-in NPT, female
Air-Cooled Condenser		Standard	Standard	-	-
Water-Cooled Condenser		Requires 2 to 4 gpm 12 °C to 22 °C	Requires 2 to 4 gpm 12 °C to 22 °C	Requires 3 to 6 gpm 12 °C to 22 °C	Requires 3 to 6 gpm 12 °C to 22 °C

Specifications are dependent on fluid choice, fluid temperatures, ambient temperature and electrical supply. Please contact SP to further discuss your requirements. RC311 and RC311GWT are water-cooled only.

¹Control Stability is process dependent

²At 50 Hz, capacity may diminish up to 17% based on conditions and application

³Water connections are for water-cooled units only



935 Mearns Road, Warminster, PA 18974 USA | scientificproducts.com +1-800.523.2327 | hello@spindustries.com

©2022 SP Industries, Inc. SP reserves the right to change specifications without notice. All trademarks referenced in this document are registered Part Nun or common law trademarks of SP Industries, Inc. or its affiliates in the U.S. and/or other countries.