

Dimensional Data

Width	35 in (89 cm)
Depth	29 in (74 cm)
Height	37 in (94 cm)
Approximate weight range	280-495 lb (127-225 kg)

Additional Information


Construction	316L stainless steel condenser chamber
Vacuum pump (required, not incl.)	Two-stage rotary vane
Defrost type	Hot gas
Refrigerant type	CFC-free

Optional Components



Vertical Manifold Assembly

4-, 8-, 12-, 24- or 48-Port Stainless Steel Vertical Manifold



Horizontal "T" Manifold

1-, 2- or 3-Tier Stainless Steel Horizontal "T" Manifold

Note: Additional standard and custom accessories are available, along with flask adapters, glassware and other components. Contact SP for more information.

* Receptacle configuration for 50 Hz units varies by country.

† The specified Maximum Ice Condensing Capacity in 24 Hours and Maximum Deposition Rate are based on the process of freeze-drying water as aggressively as possible.

The freeze dryer's ability to collect ice at an hourly rate or over a specified period will always be application dependent.

Note: The refrigerants and insulating foam contain fluorinated greenhouse gases.