



IRTA DOSIFICACIO I TECNOLOGIA, SL
Tel. +34 93 796 26 68
Carrer dels Remences, 44 P. I. Vallveric
08304 Mataró (BARCELONA)

www.i-dositecno.com



MI-R

TECHNICAL FEATURES:

Filling and capping machine, specially designed for nasal spray bottles and screw cap with different format.

Under cGMP, FDA standards and CE and UL compliance.

Máquina llenadora de envases para spray nasal y cerradora de tapones de rosca para diferentes formatos.

Fabricada bajo estándares cGMP, FDA cumpliendo con las normativas CE y UL.

TECHNICAL DATA:

Output up to: 6.000 u/h. / *Producción hasta: 6.000 ud/h.*

Filling dose Simplex: 0,1 ml- 250 ml. / *Dosificado Simplex: 0,1 ml - 250 ml.*

Documentation: DQ.

Vials rejection / *Rechazo vial.*

Signalling lights / *Señales luminosas.*

CFR21 part 11 / *CFR21 parte 11.*

OPTIONAL EQUIPMENT:

IPC 100% weight control / *IPC 100% Control de peso.*

100% torque control in capping unit./ *100% control de cierre de tapón mediante control de par.*


Laminar Flow (LF) / *Flujo Laminar (FL).*

IQ, OQ.

VPC and NVPC integration / *Integración de VPC y NVPC.*

Petri dish integration / *Integración del plato Petri.*

UPS



1- DOSI-PROX

An i-DOSITECNO trained local partner provides worldwide direct local services without any language barrier, and on the same time zone.

Un representante autorizado de i-Dositecno brinda servicios locales directos en todo el mundo sin barreras idiomáticas y en la misma zona horaria.



2- DOSI-EYE

Adopting the latest technology, i-DOSITECNO is the world's first company to provide DOSI-EYE remote services.

Every i-DOSITECNO machine is equipped with a pair of DOSI-EYE goggles, which enable our customer service engineers to provide remote video communication in real time, allowing technical guidance and troubleshooting assistance, hence solving any incident in a quick and effective manner.

Adoptando la última tecnología, i-Dositecno es la primera compañía del mundo en proporcionar el servicio remoto DOSI-EYE.

Cada máquina está equipada con gafas DOSI-EYE, lo que permite a nuestros ingenieros proporcionar al cliente comunicaciones de video remotas en tiempo real asistiendo y solucionando cualquier problema o avería de manera rápida y efectiva.



3- DOSI-GO

i-DOSITECNO commits to dispatch service engineer worldwide in less than 24 hours after receiving a service request call from our customers.

i-Dositecno se compromete a enviar un ingeniero a cualquier parte del mundo en menos de 24 horas tras recibir la petición de servicio por parte de nuestro cliente.