

miVac™ 2 Series

Modular Concentrator



DNA Concentrator

Max. rotor speed	1465 rpm
Nominal sample load g force	up to 250 g
Temperature control range	30°C to 80°C
Dimensions (W x D x H)	360 x 602 x 300 mm
With catchpot	501 x 602 x 300 mm
Weight (Approx)	41.5 kg

Duo Concentrator

Max. rotor speed	1465 rpm
Nominal sample load g force	up to 250 g
Temperature control range	30°C to 80°C
Dimensions (W x D x H)	360 x 424 x 300 mm
Weight (Approx)	21 kg

Quattro Concentrator

Max. rotor speed	1130 rpm
Nominal sample load g force	up to 250 g
Temperature control range	30°C to 80°C
Dimensions (W x D x H)	480 x 594 x 300 mm
Weight (Approx)	35 kg

Duo Pump

Ultimate vacuum [§]	<10.0 mbar (7.5 torr)
Peak performance	1.0 / 1.1 m³/h
Dimensions (W x D x H)	215 x 389 x 300 mm
Weight	19 kg

Quattro Pump

Ultimate vacuum [§]	<3.0 mbar (1.5 torr)
Peak performance	2.0 / 2.2 m³/h
Dimensions (W x D x H)	215 x 389 x 300 mm
Weight	18 kg

Super Vacuum Pump

Ultimate vacuum [§]	<0.05 mbar (7.5 torr)
Peak performance	6.2 m³/h
Dimensions (W x D x H)	282 x 432 x 302 mm
Weight	26.6 kg

Speed Trap

Type	Single-stage vapour compression
Refrigerant gas	R1270
Refrigerant charge	28 g
CO ₂ E	<0.001 tonnes
Ultimate low temperature [†]	-50 °C
Max. pressure (PS)	30 bar
Refrigeration ACC range	110 °C
Condenser jar capacity	1 litre (maximum usable volume: 800ml)
Dimensions (W x D x H)	185 x 595 x 507 mm
Weight	34 kg

Pressure Controller

Vacuum range	1-1100 mbar
Dimensions (W x D x H)	195 x 178 x 105 mm
Weight (Approx)	1.5 kg

Emissions

Noise (@ 1 meter)	74 dBA
-------------------	--------

Actual dBA value may vary depending on equipment setup / composition

Electrical Supply

100 V	50 Hz
100 V	60 Hz
120 V	60 Hz
220 V	60 Hz
230 V	50 Hz
Max. supply input	1500 A

Storage/ Transportation Environment*

Ambient Temperature [‡]	0°C to 40°C
Relative humidity	10-80% non-condensing

**Store upright at all times.*

Operational Environment*

Ambient temperature	150°C to 30°C
Relative humidity	10-80% non-condensing
Altitude	Sea-level to 1600 m
Min. ventilation air-gap	50 mm
Installation environment	Indoor only

Static- dissipative laboratory or similar

[§] Ultimate low temperature; operational values may vary

[†] Ultimate vacuum; operational values may vary

[‡] -10°C permissible during transport (only)



935 Mearns Road, Warminster, PA 18974 USA | scientificproducts.com
+1-800.523.2327 | hello@scientificproducts.com

©2023 SP Industries, Inc. SP reserves the right to change specifications without notice. All trademarks referenced in this document are registered or common law trademarks of SP Industries, Inc. or its affiliates in the U.S. and/or other countries. *Other specific brand or product names may be registered trademarks of their respective companies.

Rev 001, 01/25