

## **Laboratory Glassware Washers**

# **Specifications**

### **Dimensional Data (Nominal)**

## Undercounter

Width	610 mm (24.00 in)
Depth	580 mm (22.88 in)
Height	875 to 924 mm (34.44 to 36.38 in)
Weight	53 kg (116.80 lbs)
Volume	170 L (6 ft <sup>3</sup> )

#### Freestanding

Width	610 mm (24.00 in)
Depth	616 mm (24.25 in)
Height	978 mm (38.50 in)
Weight	62 kg (136.50 lbs)
Volume	170 L (6 ft <sup>3</sup> )

#### Vertical SpaceSaver™

The state of the s		
Width	670 mm (26.38 in)	
Depth	592 mm (23.31 in)	
Height	1994 mm (78.50 in)	
Weight	168 kg (370.00 lbs)	
Volume	340 L (12 ft <sup>3</sup> )	

## Water Consumption per Chamber Minimum (excluding Rinse Only cycle)

Тар	2.4 Gal (9 L)
DI	0.8 Gal (3 L)

## Maximum (excluding Rinse Only cycle)

Тар	4 Ga	al (15 L)	)
DI	1.6	Gal (6 L	)



#### **Material of Construction**

Chamber Material	304 Stainless Steel
Exterior Material	304 Stainless Steel
Certification	UL61010-1 & CSA 22.2 61010-1
Electrical	120V, 60Hz, 15A, with NEMA 5-15P power cord

## Warranty

SP Hotpack washers and accessories are covered for a period of 12 months from date of shipment. For full details contact your SP Hotpack area representative.

#### Note:


The SP Hotpack recommendation is for at least 673 mm (26.5 in) of clearance in front of the unit for the door extension and at least 51 mm (2 in) on each side of the door.





## **Laboratory Glassware Washers**

### **Undercounter Glassware Washers**



## **Vertical Spacesaver Glassware Washers**

